

CB1 - Premium Rounded Cabinets

The CB1 Cabinet comes in both large and small variants, designed to house all AED models from major manufacturers while maintaining the smallest wall footprint possible. Strobed and non alarmed versions are available for both variants of cabinets.



Our CB1 line of cabinets stand alone in the AED industry with a unique set of premium innovations combining a sleek aesthetic, thoughtful safety features, structural integrity and AED compatibility with no compromise. Its stand-out features are not limited to its trademark rounded safety corners, but also include a recessed handle, concealed hinges, magnetic door seal, 1.1mm cold rolled steel fabrication, and an impressive textured powder coating finish. The 80-120db keyed alarm system housed in a durable internal enclosure offers security for the AED.

In addition to these impressive product features, our cabinets are designed to stand the test of time; with a standard lifetime warranty, and the best packaging in the industry, our cabinets are made for life.

Product Photos



Product Dimensions

CB1-L Premium Large

- Length: 16"
- Width: 8-3/8"
- Height: 14-5/8"

CB1-S Premium Small

- Length: 14-3/4"
- Width: 6-3/4"
- Height: 11-5/8"

Engineered For Safety

Our CB1 Premium cabinets feature fully-welded rounded corners, providing a safe and aesthetically pleasing storage solution for your AED and First Responder Kits. Ideal for high traffic areas, especially those where young children are present.

